INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only
CE 0086

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name: **Azure[®] Anterior Cervical Plate System**

 Click directory below for desired language

English 	2-4
Español 	5-7

INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only
CE 0086

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

Azure[®] Anterior Cervical Plate System

Description:

The Azure Anterior Cervical Plate System consists of an assortment of nonsterile, single use, titanium alloy (Ti-6Al-4V ELI per ASTM F136) with nitinol (per ASTM F2063) plates and screws that allow a surgeon to build a temporary anterior cervical implant construct. The plate is attached to the anterior aspect of the vertebral body by means of screws to the cervical spine. The system includes the necessary instrumentation to assist in the surgical implantation of the devices.

Indications for Use:

The Azure Anterior Cervical Plate System is a temporary implant intended for anterior fixation to the cervical spine from C2 to C7 and indicated for:

1. Degenerative disc disease (defined as back pain of discogenic origin with degenerative disc confirmed by patient history and radiographic studies).
2. Spondylolisthesis.
3. Trauma (i.e., fracture or dislocation).
4. Spinal stenosis.
5. Deformities (i.e., scoliosis, kyphosis, and/or lordosis).
6. Tumor.
7. Pseudoarthrosis.
8. Revision of previous surgery.

Contraindications:

The Azure Anterior Cervical Plate System is contraindicated in patients with a systemic infection, with a local inflammation at the bone site or with rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis. Do not use this system in patients with known or suspected metal allergies. Use of the system is also contraindicated in patients with any other medical, surgical or psychological conditions that would preclude potential benefits of internal fixation surgery such as the presence of tumors, congenital abnormalities, elevation of sedimentation rate unexplained by other disease, elevation of white blood cells or a marked shift in white blood cell differential count.

Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Early or late loosening of any or all of the components.
2. Disassembly, bending, and/or breakage of any or all of the components.
3. Foreign body (allergic) reaction to implants, debris, corrosion products, and graft material, including metallosis, straining, tumor formation, and/or auto-immune disease.
4. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain.
5. Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
6. Infection.
7. Vertebral body fracture at, above, or below the level of surgery.
8. Loss of neurological function, including paralysis (complete or incomplete).
9. Non-union, delayed union.
10. Pain, discomfort, or abnormal sensations due to the presence of the device.
11. Hemorrhage.
12. Cessation of any potential growth of the operated portion of the spine.
13. Death.

Note: Additional surgery may be necessary to correct some of these anticipated adverse events.

Warnings and Precautions:

The surgeon should be aware of the following when using this implant system:

1. The Azure Anterior Cervical Plate System is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical spine.
2. To optimize bony union, perform an anterior microdiscectomy or corpectomy as indicated.
3. To facilitate fusion, a sufficient quantity of autograft or allograft material should be used.
4. Always orient the plate along the midline of the spine.
5. Excessive torque applied to the screws when seating the plate may strip the threads in the bone.
6. Do not leave screws partially lagged (directly below the plate locking mechanism) or in contact with the locking mechanism. Screws left proud (partially lagged) may result in an upward force on the locking mechanism when tightening adjacent screws which may cause failure of the plate locking mechanism. Ensure screws are fully lagged or left above the locking mechanism before final tightening of the entire construct.

7. Failure to place screws within the intended angle range may result in plate locking mechanism fracture or screw back-out past the plate locking mechanism.
8. Failure to achieve arthrodesis will result in eventual loosening and failure of the device construct.
9. Do not reuse implants; discard used, damaged, or otherwise suspect implants.
10. Implants are single use only. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use implants that come in contact with body fluids.
11. The health care provider must thoroughly clean the instruments prior to steam sterilization.
12. NON-STERILE; the plates, bone screws and instruments are non-sterile, and therefore, must be sterilized before each use.
13. Do not combine Azure Anterior Cervical Plate System implantable components with those from any other system or manufacturer.
14. Do not combine implants of dissimilar metals, as it may result in galvanic corrosion.

MRI Compatibility Information:

The Azure Anterior Cervical Plate System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The system has not been tested for heating or migration in the MR environment.

Cleaning:

Implants are provided clean, but not sterile. Once an implant comes in contact with any human tissue or bodily fluid it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned and sterilized prior to each use. Reprocess all instruments as soon as is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

From Point of Use:

Wherever possible, do not allow blood, debris, or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

Preparation for Cleaning:

1. There are two instruments in the Azure system, the Plate Holder (86-9003) and the DTS Guide (86-9006), which require disassembly prior to cleaning. No other instruments within the Azure system require disassembly prior to cleaning.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to manual or automated cleaning process.
3. Any instruments with moving parts (i.e., knobs, triggers, hinges) should be separated and activated to open position to allow better access of the cleaning fluid to the difficult to clean areas.
4. Use a soft cloth or plastic bristle brush to remove any visible soil from the instruments.

Disassembly Instructions for the Plate Holder (86-9003):

Remove the inner shaft from the plate holder shaft by pulling axially until the inner shaft is free.

Disassembly Instructions for the DTS Guide (86-9006):

Remove the inner shaft from the guide shaft by pulling axially until the inner shaft is free.

Manual Cleaning:

1. After the disassembly of the instruments, use a soft cloth and/or a soft plastic bristle brush to remove any visible soil from the outside and inside of the instruments.
2. Prepare Vesphene[®] Ise agent at the use-dilution recommended by label direction (1 ounce per gallon) or 1 mL of Vesphene[®] Ise to 128 mL of potable tap water per manufacturer recommendations.
3. Bathe instruments in prepared room temperature solution as recommended by the detergent manufacturer.
4. Manually agitate instruments in Vesphene[®] solution for 15 minutes.
5. Scrub instruments with soft plastic bristle brush if visible soil is noted and use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil.
6. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
7. Hang dry the device.
8. Visually inspect the instruments for visible soil.
9. If visible soil is noted, repeat the steps listed above.

Automated Cleaning:

1. After disassembly of the instruments use a soft cloth or soft plastic bristle brush to remove any visible soil from the outside of and inside of the instruments. Use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil from the inner lumens prior to automated cleaning.
2. Set up the washer/disinfector detergent dose as ½ ounce of Endozime AW Plus® per gallon of water or according to manufacturer recommendations.
3. Remove the instruments from the instrument case and place them into the washer baskets.
4. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
 - a. Pre-Wash 1: cold potable water, 2 minutes
 - b. Enzyme/Detergent treatment
 1. Spray, 20 seconds
 2. Soak, 1 minute
 3. Rinse cold potable water, 15 seconds
 4. Rinse cold potable water, 15 seconds
 - c. Wash ≥ 65°C, 2 minutes using a Endozime AW Plus®
 - d. Rinse 1: hot potable water, 15 seconds
 - e. Rinse 2: hot potable water, 15 seconds
 - f. Rinse 3: hot potable water, 15 seconds
 - g. Rinse 4: hot potable water, 15 seconds
 - h. Thermal rinse ≥ 93°C for 1 minute
 - i. Heated USP Purified Water Rinse 1: re-circulating 10 seconds
 - j. Heated USP Purified Water Rinse 2: non re-circulating 10 seconds
 - k. Drying: 7 minutes, 115°C
6. Visually inspect the instruments for visible soil.
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach and/or other alkaline cleaners may cause damage, particularly to instruments; these solutions should not be used.

Note: Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix if instruments are broken, discolored, corroded, have cracked components, pits, gouges, or otherwise found defective. Do not use defective instruments.

Assembly/Inspection Instructions for the Plate Holder (86-9003):

1. Align both the plate holder shaft and the inner shaft.
2. Insert the inner shaft into the plate holder shaft until it clicks once.
3. After assembly, ensure the shaft slides smoothly between the locked and unlocked positions.
4. If the shaft does not slide smoothly between the locked and unlocked positions, disassemble the instrument and check for debris.
5. If debris is found, repeat cleaning and attempt assembly instructions again.

Assembly/Inspection Instructions for the DTS Guide (86-9006):

1. Align both the DTS guide shaft and the inner shaft.
2. Insert the inner shaft into the guide shaft until it clicks once.
3. After assembly, ensure the shaft slides smoothly between the locked and unlocked positions.
4. If the shaft does not slide smoothly between the locked and unlocked positions, disassemble the instrument and check for debris.
5. If debris is found, repeat cleaning and attempt assembly instructions again.

Sterilization:

The Azure Anterior Cervical Plate System implants and instruments are supplied NON-STERILE. Prior to use, all implants and instruments should be placed in the instrumentation/implant case which will be wrapped in an FDA cleared sterilization wrap or placed in a rigid sterilization container and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Sterilization in Orthofix Instrument Cases with Blue Wrap:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Preconditioning: 4 Pulses
Drying time: 30 minutes		Exposure time: 4 minutes
Double wrapped		Drying time: 30 minutes
(FDA cleared wrap recommended)		Double wrapped
		(FDA cleared wrap recommended)

Sterilization in Rigid Sterilization Containers:

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer instructions.

Select the appropriate rigid sterilization container (either filtered bottom or solid bottom) to properly enclose the Orthofix instrumentation case(s) (recommended 23¼" long x 11¼" wide container). The following sterilization cycle was validated:

Method: Steam
Cycle: Prevac
Temperature: 270°F (132°C)
Preconditioning: 4 Pulses
Exposure time: 4 minutes
Drying time: 30 minutes

Validation and routine monitoring should be performed as per ANSI/AAMI ST79 Comprehensive guide to steam sterilization and sterility assurance in health care facilities. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Plate Holder (86-9003) Use Instructions:

1. Once the plate holder shaft is attached to the plate holder square in the plate, push the inner shaft to lock.
2. Pull the inner shaft to unlock.

DTS Guide (86-9006) Use Instructions:

1. Once the guide shaft is attached to the plate holder square in the plate, push the inner shaft to lock.
2. Pull the inner shaft to unlock.

Physician Information:

Patient Selection:

Patient selection is an extremely important factor in the success of implant procedures. It is important that the candidates are carefully screened and the optimal therapy selected.

Preoperative:

1. Carefully screen the patient, choosing only those that fit the indications described above.
2. Care should be exercised in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Store away from corrosive environments.
3. An adequate inventory of implant sizes should be available at the time of surgery.
4. All components and instruments should be cleaned and sterilized prior to each use. Additional sterile components should be available in case of an unexpected need.

Intraoperative:

1. Instructions should be carefully followed.
2. Extreme caution should be used around the spinal cord and nerve roots.
3. The implant surface should not be scratched or notched since such actions may reduce the functional strength of the construct.
4. Bone grafts must be placed in the area to be fused such that the graft fits snugly against the upper and lower vertebral bodies.
5. Bone cement should not be used as it will make removal of the components difficult or impossible.
6. Before closing soft tissue, check each screw to make sure that none have loosened.

Postoperative:

1. Detailed instructions should be given to the patient regarding care and limitations, if any.
2. To achieve maximum results, the patient should not be exposed to excessive mechanical vibrations. The patient should not smoke or consume alcohol during the healing process.
3. The patient should be advised of their limitations and taught to compensate for this permanent physical restriction in body motion.
4. If a non-union develops or if the components loosen, the devices should be revised or removed before serious injury occurs. Failure to immobilize the non-union or a delay in such will result in excessive and repeated stresses on the implant. It is important that immobilization of the spinal segment be maintained until fusion has occurred.
5. The implants are temporary internal fixation devices. Internal fixation devices are designed to stabilize the spine during the normal healing process. After the spine is fused, the devices serve no functional purpose and should be removed.
6. In most cases removal is indicated because the implants are not intended to transfer or support forces developed during normal activities.
7. If the device is not removed following completion of its intended use, complication may occur as follows:
 - a. Corrosion, with localized tissue reaction or pain.
 - b. Migration of implant position resulting in injury.
 - c. Risk of injury from postoperative trauma.
 - d. Bending, loosening and/or breakage, which could make removal impractical or difficult.
 - e. Pain, discomfort or abnormal sensations due to the presence of the device.
 - f. Possible increased risk of infection.
 - g. Bone loss caused by stress shielding.

Patient Information:

The temporary internal fixation devices used in your recent spinal surgery are metallic implants that attach to the bone and aid in the healing of bone grafts. These implants have been shown to be valuable aids to surgeons in the treatment of bony fusions. These devices do not have the capabilities of living bone. Intact living bone is self-repairing, flexible and occasionally breaks and/or degrades. The anatomy of the human body places a size limitation on any artificial fixation device used in surgery. The maximum size limitation increases the chances of the mechanical complications of loosening, bending, or breaking of the devices. Any of these complications could result in the need for additional surgery. Accordingly, it is very important that you follow the recommendations of your physician. Use braces as instructed. By following these instructions, you can increase your chances of a successfully result and reduce your risk of injury and/or additional surgery.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The Azure Anterior Cervical Plate System instruments and implants are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700, or by e-mail at complaints@orthofix.com.

Further information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The implants, instruments and/or packaging material for the Azure Anterior Cervical Plate System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Rx Only	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use	 Manufacturer
	Orthofix.com/IFU	 Authorized Representative
	Single Use Only Do Not Reuse	 Serial Number
	Catalogue Number	 Lot Number
	Provided Non-Sterile	

Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

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