INSTRUCTIONS FOR USE
Important Information – Please Read Prior to Use

Rx Only
CE

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
QUADRx[®] Lateral Retractor System

Click directory below for desired language

English **EN** **2-3**

Español **ES** **4-5**

INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only
CE

Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

QUADRx® Lateral Retractor System

Description:

The QUADRx Lateral Retractor System consists of manual surgical instruments intended to aid the surgeon's access and visualization of the surgical area during spinal procedures.

The QUADRx Lateral Retractor System includes stainless steel detachable blades in multiple lengths, fiber optic lighting, and table attachments to connect to the side rail of the operating room table.

Indications for Use:

The QUADRx Lateral Retractor System consists of manual surgical instruments intended to aid the surgeon's access and visualization of the surgical area during spinal procedures.

Contraindications:

The QUADRx Lateral Retractor System is contraindicated for use in patients:

1. Morbid obesity.
2. Mental illness.
3. Alcoholism or drug abuse.
4. Pregnancy.
5. Metal sensitivity/allergies.
6. Severe osteopenia.
7. Patients unwilling or unable to follow post-operative care instructions.
8. Any circumstances not listed under the heading Indications for Use.

Potential Adverse Events:

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Device component fracture.
2. Neurological injury.
3. Vascular or visceral injury.
4. Foreign body (allergic) reaction to instruments, debris, corrosion products, including metallosis, straining, tumor formation, and/or autoimmune disease.
5. Infection.
6. Hemorrhage.
7. Cessation of any potential growth of the operated portion of the spine.
8. Death.

Note: Potential risks identified with the use of the device system may require additional surgery.

Warnings and Precautions:

1. The QUADRx Lateral Retractor System is sold non-sterile and therefore must be sterilized before use.
2. Care should be exercised in the handling and storage of instruments. Instruments should not be scratched, notched or otherwise damaged since such actions may reduce functional performance. Store away from corrosive environments.
3. The retractor should be assembled prior to surgery. An adequate inventory of blades should be available during surgery in case a size other than those expected to be used is needed.
4. All components and instruments should be cleaned and sterilized prior to use. Additional sterile components should be available in case of unexpected need.
5. If used around the spinal cord and nerve roots, extreme caution should be taken.
6. The reusable fiber optic light cables are designed for use with 300 watt xenon illuminators using the provided ACMI adapter cable. Do not use light sources rated higher than 300 watts or any cables other than the provided ACMI adapter cable. Use of higher watt sources or cables other than the provided ACMI adapter cable could result in overheating, causing product failure and patient injury.
7. Do not operate the light source and adapter cable without the reusable fiber optic light cable attached. Without the reusable fiber optic light cable, the output from the adapter cable is extremely bright, hot, and may cause burns, ignite drapes/gowns, or temporarily blind vision.
8. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use instruments that come in contact with body fluids.

Cleaning:

All instruments must be thoroughly cleaned prior to use and after each re-use of the instruments. Reprocess all instruments as soon as it is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

None of the instruments in the system require disassembly prior to cleaning.

From Point of Use:

Wherever possible, do not allow blood, debris, or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

Preparation for Cleaning:

1. Any instruments with moving parts (e.g., knobs, triggers and hinges) should be placed in the open position to allow better access of the cleaning fluid to the difficult to clean areas.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process.
3. Use a soft cloth or plastic bristle brush to remove any visible soil from the instruments.

Manual Cleaning:

1. Because none of the instruments can be disassembled, use a soft cloth and/or a soft plastic bristle brush to remove any visible soil from the outside and inside of the instrument.
2. Prepare Vesphene® Ilse agent at the use-dilution recommended by label direction – 1 ounce per gallon or 1 mL of Vesphene® Ilse to each 128 mL of potable tap water per manufacturer's recommendations.
3. Bathe instruments in prepared room temperature solution as recommended by the detergent manufacturer.
4. Manually agitate instruments in Vesphene® Ilse solution for 15 minutes.
5. Scrub instruments with a soft plastic bristle brush if visible soil is noted and use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil.
6. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
7. Hang dry the device.
8. Visually inspect the instruments for visible soil.
9. If visible soil is noted, repeat the steps listed above.

Automated Cleaning:

1. Use a soft cloth or soft plastic bristle brush to remove any visible soil from the outside and inside of the instruments. Use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil from the inner lumens prior to automated cleaning.
2. Set up the washer/disinfectant detergent dose as ½ ounce of Endozime AW Plus® per gallon of water or according to the manufacturer recommendations.
3. Remove the instruments from the instrument case and place them into the washer baskets.
4. Orient instruments in the automated washer's carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
 - a. Pre-Wash 1: cold potable water, 2 minutes
 - b. Enzyme/Detergent treatment:
 1. Spray, 20 seconds
 2. Soak, 1 minute
 3. Rinse cold potable water, 15 seconds
 4. Rinse cold potable water, 15 seconds
 - c. Wash ≥ 65°C, 2 minutes using Endozime AW Plus®
 - d. Rinse 1: hot potable water, 15 seconds
 - e. Rinse 2: hot potable water, 15 seconds
 - f. Rinse 3: hot potable water, 15 seconds
 - g. Rinse 4: hot potable water, 15 seconds
 - h. Thermal rinse ≥ 93°C, 1 minute
 - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
 - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
 - k. Dry at 115°C, 7 minutes
6. Visually inspect the instruments for visible soil
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

Note: Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix if instruments are broken, discolored, corroded, have cracked components, pits, gouges, or otherwise found defective. Do not use defective instruments.

Sterilization:

QUADRx Lateral Retractor System instruments are supplied NON-STERILE. Prior to use, all instruments should be placed in the instrument case which will be either wrapped in an FDA cleared sterilization wrap or placed in a rigid sterilization container and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Sterilization in Orthofix Cases with Blue Wrap:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Exposure time: 4 minutes
Drying time: 30 minutes		Drying time: 30 minutes
Double wrapped		Double wrapped
(FDA cleared wrap recommended)		(FDA cleared wrap recommended)

Sterilization in Rigid Sterilization Containers:

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer instructions.

Select the appropriate rigid sterilization container with either a filtered or solid bottom to properly enclose the Orthofix case (recommended 23¼" long x 11¼" wide container). The following sterilization cycle has been validated:

Method: Steam
 Cycle: Prevac
 Temperature: 270°F (132°C)
 Exposure time: 4 minutes
 Drying time: 30 minutes

Note: Rigid sterilization containers with solid bottoms cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The QUADRx Lateral Retractor System instruments are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at complaints@orthofix.com.


Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The instruments and/or packaging material for the QUADRx Lateral Retractor System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Rx Only		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use		Catalogue Number
	Orthofix.com/IFU		Manufacturer
	Single Use Only Do Not Reuse		Authorized Representative
	Provided Non-Sterile		Lot Number

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